

Open-Top Polyclear Tubes

Tubes Size	PN Box of 50	Tubes Capacity ml	For Beckman Rotor	For Sorvall/Kendro Rotor	For Kontron Rotor	Equivalent for Beckman	Equivalent for Sorvall/Kendro
3/16x1-15/16" (5x41mm)	7009	.80 nom .67act	SW60Ti (4) SW65Ti, SW55Ti SW50.1, SW50	TH-660 (4) AH-650	---	344090	---
7/16x1-3/8" (11x34mm)	7011	2.2	TLS55	RP55-S	---	347356	45315 45242
7/16x2-3/8" (11x60mm)	7010	4.4	SW60Ti SW56Ti	TH-660	TST-60.4 TFT-80.4	344062	---
1/2x1-1/4" (13x32mm)	7020	3.0	90Ti, 80Ti, 75Ti, 70.1Ti, 50Ti, 65, 50 & 40 Adapter required	---	---	344092	---
1/2x 1-5/8" (13x41mm)	7021	4.0	90Ti, 80Ti, 75Ti, 70.1Ti, 1Ti 50Ti, 65, 50 & 40 Adapter required	---	---	344093	---
1/2x2" (13x51mm)	7022	5.0	SW65Ti, SW55Ti SW50.1, SW50	AH-650	TST55.5	344057	---
1/2x2-1/2" (13x64mm)	7023	6.5	50.3Ti, 40.3 40.2, 50.4T1 (1)	---	TFT45.6	344088	---
9/16x3-1/2" (14x89mm)	7030	12.0	SW41Ti	TH-641	---	344059	06752
9/16x3-3/4" (14x95mm)	7031	12.9	SW40Ti	---	TST41.14	344060	---
5/8x2" (16x51mm)	7043	8.0	SW30.1	---	---	346600	---
5/8x2-1/2" (16x64mm)	7040	10.2	Ty50	---	---	344089	---
5/8x3" (16x76mm)	7041	12.6	90Ti, 80Ti, 70.1Ti 75Ti, 65, 50Ti 40, 50(2)	T-1270, T-875 T-865, A-1256	TFT80.13, 75.13 65.13, 60.13 50.13, 32.13	344085	---
5/8x4" (16x102mm)	7042	17.0	SW28.1Ti SW27.1Ti SW32.1	AH-629/17 Surespin 630/17	TST28.17	344061	06750
1x2" (25x51mm)	7053	20.0	SW30	AH-629/20	---	346598	---
1x3" (25x76mm)	7050	32.5	SW25.1	---	---	344063	---
1x3-1/2" (25x89mm)	7052	38.5	SW32, SW27, SW28 70Ti, 60Ti, 55.2Ti 50.2Ti, 42.1, 30	AH-629/36, T1250 Surespin 630/36	TFT70.38, 65.38, 50.38, 55.38 TST28.38	344958	75000471
1-1/2x4" (38x102mm)	7060 (3)	94.0	45Ti, 35, 21	---	TFT45.94	345777	---

(1) May be used with adapters in the following rotors: Types 90Ti, 80Ti, 75Ti, 70.1Ti, 65, 60Ti, 55.2Ti, 50.2Ti, 50Ti, 45Ti, 42.1, 40, 35 and 21.

(2) May be used with adapters in the following rotors: Types 70Ti, 60Ti, 55.2Ti, 50.2Ti, 45Ti, 42.1, 35, 30 and 21.

(3) PN 7060 is a box of 25 tubes.

(4) In SW60Ti and TH-660 rotors (11mm dia) use Adapter PN 4109 (Set of 6 adapters);

in SW65Ti, SW55Ti, SW50.1, SW50 rotors (13mm dia) use Seton Adapter PN 4114 (Set of 6).